

## In gel protein digestion using trypsin as protease, and extraction of liberated peptides

### In gel digestion (day 1)

#### Silver removal (only silver stained gels):

Cut gel bands/spots into 1 mm cubes in eppendorf tubes. Wash once with 500  $\mu$ l water (incubate at room temperature (RT) for 20 min. in an 1.5 ml Eppendorf Thermomixer comfort, Eppendorf AG, Germany). Add 200  $\mu$ l destain solution ([see right panel](#)), and incubate at RT in Eppendorf mixer for 5 min (slow agitation).

Remove supernatant, and wash gel pieces by adding 100  $\mu$ l MilliQ water. Incubate at RT in Eppendorf mixer for 5 min (slow agitation), and discard supernatant. Repeat 4 times.

#### Destain solution:

Mix equal volumes of 30 mM potassium ferricyanide (10 mg  $K_3Fe(CN)_6$ /ml MilliQ water) and 100 mM sodium thiosulfate (16 mg  $Na_2S_2O_3$ /ml MilliQ).

#### Gel washing:

Cut gel bands/spots into 1mm cubes in eppendorf tubes. Add 50-100  $\mu$ l wash solution ([see right panel](#)) to each sample, and incubate at RT for 20 min. in the Eppendorf mixer. Remove supernatant, and repeat wash once.

Discard supernatant, dry gel pieces preferably in vacuum, a "Rotavapor" (or by adding 50  $\mu$ l ACN and shake for 2 min. Remove ACN. The gel pieces should now be white and sticky).

#### Wash solution:

Add 250  $\mu$ l 1M Ambic (frozen stock solution of 1M ammoniumbicarbonate in MilliQ water, 80 mg/ml) to 4750  $\mu$ l MilliQ water and 5 ml ACN (acetonitrile, HPLC grade)

#### In gel reduction og alkylation (Cys):

Reduce cysteins by adding 50  $\mu$ l 10 mM DTT (DiThioTreitol from Amersham Biosciences, #171318-02) to the dried gel pieces ([see right panel](#)), and incubate at 56 °C for 45 minutes.

Cool samples, and remove DTT solution. Immediately add 50  $\mu$ l 55 mM IAA (iodoacetamide, Sigma Aldrich, I-6125) for cystein alkylation ([see right panel](#)), and incubate in the dark at room temperature for 30 min. Remove IAA solution, wash twice as described above, and dry gel pieces in vacuum, a "Rotavapor".

#### 10 mM DTT in 100 mM Ambic:

Add 10  $\mu$ l 1M DTT (Frozen stock solution of 154 mg DTT/ml MilliQ water) to 890  $\mu$ l MilliQ water and 100  $\mu$ l 1M Ambic.

#### 55 mM IAA in 100 mM Ambic:

Add 10 mg IAA to 900  $\mu$ l MilliQ water and 100  $\mu$ l 1M Ambic.



In gel protein digestion;

*Coomassie spots:* Add 20 - 40µl 6ng/µl Trypsin Porcine (from Promega, #V 511A) to each sample (*see right panel*), and rehydrate on ice for 30 min.

*Silver and Sypro spots:* Add 20 - 40µl 3ng/µl Trypsin Porsine (from Promega, #V 511A) to each sample (*see right panel, but using 5µl Trypsin stock solution*), and rehydrate on ice for 30 min.

Incubate samples for 16 hours at 37 °C in a hot cabinet.

Digestion buffer:

Add 50 µl 1M Ambic, 10µl 0.1M CaCl<sub>2</sub> and 50 µl ACN, to 890 µl MilliQ water.

6ng/µl Trypsin:

Mix 10 µl Trypsin Promega Porsine stock solution (100ng/µl dissolved in 50 mM acetic acid) and 160 µl digestion buffer (*see above*).

**Peptide extraction (day 2)**

Extraction of peptides:

Cool and spin samples. Pull off and save supernatant in **new eppendorf tube**.

Add 30-50 µl 5% FA (formic acid), and incubate at room temperature for 20 min in the eppendorf mixer. Pull off supernatant and pool with the first extraction.

Add 30-50 µl 60% ACN/0.1 % FA to gel samples, and incubate for 20 min. in the eppendorf mixer. Pull off supernatant and pool with the two former extractions.

Vacuum dry solution in a Rotavapor (Concentrator 5301 from Eppendorf AG, Hamburg, Germany) till it remains 10-15 µl samples. If necessary, add 15 µl 0.1% TFA and vacuum dry to about the same volume as above (samples should not contain ACN).