


Filtersets available on the Nikon Eclipse TE2000-E

Excitation: 426-446 **Emission: 460-500**



CFP - 31044v2

CFP, Chromomycin A


Excitation: 450-490 **Emission: 500 Lp**



EGFP - 41018

Calcein, Cy2, DiO, EGFP, FITC

Excitation: 490-510 **Emission: 520 Lp**




YFP - 41029

YFP, Alexa 488, DiO, FITC

Filtersets available on the Nikon Eclipse TE2000-E


Excitation: 340-380 Emission: 435-485



Dapi - 31000

Dapi, Hoechst, Calcein blue

Excitation: 465-495 Emission: 515-555



FITC - 31001

Alexa 488, Calcein, Cy2, DiO, FITC, YoYo, Yo-Pro, Eosin, YFP

Excitation: 528-552 Emission: 577-632



TRITC - 31002

Alexa 546/555, Calcium Orange, Cy3, DiI, Propidium Orange, Rhodamine/TRITC, Ethidium Bromide

Filtersets available on the Nikon Eclipse TE2000-E

Excitation: 590-650 Emission: 662-738



CY5- 41008